EN Dr. Brown's® Zero-Resistance™ Specialty Feeding System Includes the Infant-Paced Feeding Valve

Dr. Brown's® Zero-Resistance™ Specialty Feeding System was designed in collaboration with medical professionals seeking an oral feeding solution for infants identified with complex feeding challenges such as those often observed in infants with cleft lip and palate, ankyloglossia, high-arched palate, oro-neuromotor dysfunction, common or rare syndromic sequences and/or craniofacial anomalies.

The *Zero-Resistance™* Specialty Feeding System assists infants who are unable to generate adequate suction to retrieve fluid from a bottle. This bottle system allows the infant to retrieve fluid by compressing the nipple/teat. The Infant-Paced Feeding Valve (blue one-way valve) can be inserted into any Dr. Brown's narrow nipple/teat to create a "compression" nipple/teat.

Benefits of the Zero-Resistance™ Specialty Feeding System:

- Consistent and reliable nipple/teat flow rates
- Infant directed feeding vs. manually expressed boluses by feeder
- Requires minimal assembly
- Helps decrease gassiness/air intake and burping with non-squeeze method
- Preserves Vitamins A, C, and E in breastmilk.
- Easy to clean and safe for dishwasher, sanitizer, and autoclave*
- Fully reusable system
- The full bottle system and nipples/teats are BPA free and do not contain natural rubber latex.
- Exhibits a typical feeding bottle design

***Autoclave (Steam Sterilization) is not to exceed 30 cycles. Handi-Craft Company has not tested the Infant-Paced Feeding Valve (blue one-way valve) for steam sterilization/autoclaving.**

NOTE: With the Infant-Paced Feeding Valve (blue one-way valve) in place, the nipple/teat will fill with fluid as infant elevates tongue to the underside of the nipple/teat during feeding. No manual priming of the nipple/teat is necessary.

Technical Notes:

NOTE: The Zero-Resistance™ Specialty Feeding System WILL NOT FUNCTION as intended without the INSERT, RESERVOIR (Internal Vent System) and INFANT-PACED FEEDING VALVE.

- The *Zero-Resistance™* Specialty Feeding System is composed of either a 4 oz/120 mL or an 8 oz/250 mL bottle system with the internal vent system, a silicone nipple/teat and the Infant-Paced Feeding Valve.
- Eight (8) different nipple/teat flow rates are available to custom-fit all infants' needs.
- This bottle system is not intended to be a disposable system. Prior to first use, wash with soap and water. In hospital/clinical setting, consider sanitizing. In home environment, boil for 5 minutes, rinse well, and air dry. Rinse immediately after use and clean with soap and water, rinse well, and air dry before reusing. In hospital/clinical environment, consider sanitizing once every 24 hours.
- Do not warm the liquid beyond body temperature of 98.6°F (37°C).
- If you are having difficulty using the *Zero-Resistance™* Specialty Feeding System, call **Handi-Craft's Customer Service at 1-800-778-9001** or contact us at: **help@handi-craft.com**. For purchases made outside the United States and Canada, please contact your local distributor.
- For more information about the *Zero-Resistance™* Specialty Feeding System visit: **drbrownsmmedical.com/products/specialty-feeding**

Before first use: Use clean hands and clean surfaces when handling. Disassemble and clean the product and then place the components in boiling water for 5 minutes. This is to ensure hygiene. Wash all parts prior to use. Be sure the nipple/teat hole is clean.

Before each use: Pull on the nipple/teat in all directions to ensure it does not tear. Inspect carefully, especially when the child develops teeth, and replace the nipple/teat immediately if any signs of wear or damage are visible, such as cracking or change of shape. Small pieces can accidentally be inhaled and cause serious injury. Clean before each use. The nipple/teat must be clean and free of obstructions.

Cleaning: Fully disassemble the bottle system before cleaning. Be sure to remove wire cleaning brush from bottle before cleaning. Wash all parts thoroughly in hot water and detergent. The included wire brush can be used to clean the opening of the internal vent system and the Infant-Paced Valve (blue one-way valve). Rinse parts thoroughly, ensuring the nipple/teat hole is clean and clear. Air dry.

CAUTION: Do not use the wire cleaning brush to clean the nipple/teat as this will cause tears. All bottle parts, including the vent system, are dishwasher safe (top rack only), and safe for use in electric sterilizer, microwave sterilizer or boiling water.

STERILIZING BOTTLES AND NIPPLES/TEATS: You can sterilize bottle-feeding equipment by boiling, soaking in chemical/sterilizing solution or using an electric or microwave sterilizer. Use chemical/sterilizing solutions or sterilizers strictly according to the manufacturer's instructions. We recommend using the Dr. Brown's® Microwave or Electric Sterilizer. Steam can cause severe skin burns, use caution when boiling on a stove or using any sterilizing equipment. Wash your hands thoroughly with soap and water before handling sterilized items. Avoid unnecessary handling of sterilized items and equipment. Do not sterilize or microwave the wire cleaning brush.

For your child's safety and health

WARNING!

- Continuous and prolonged sucking of fluids will cause tooth decay.
- Always check food temperature before feeding.
- Throw away at the first signs of damage or weakness.
- Keep components not in use out of the reach of children.
- Never attach to cords, ribbons, laces or loose parts of clothing. The child may be strangled.
- Never use feeding nipples/teats as a pacifier/soother.
- Remove from your baby's mouth when the baby is not feeding.
- Always use this product with adult supervision.
- Glass containers may break.
- Always check the nipples/teats for signs of wear and replace when damaged.
- Do not use a microwave oven to heat the bottle as microwave ovens can create hotspots in milk/formula and/or bottle components.
- The bottle components are not toys and should not be given to children to play with or suck on.
- Do not leave the product in direct sunlight, heat or in disinfectant (sterilizing solution) for longer than recommended as this may weaken the nipple/teat.
- Replace nipples/teats after approximately 2-3 months of use.
- When used and sterilized intensively, nipples/teats may deteriorate after six weeks.
- Do not allow a child to walk or run while using any bottle.

CAUTIONS: Stop feeding before your baby has completely emptied the nipple/teat to keep your baby from swallowing air that may cause colic or other digestive problems. Be sure to contact your physician should you have any feeding problems or questions.

GLASS WARNING: Glass bottles may break. Bottles should be inspected before each use for sharp edges, fragments or flakes of glass that may have become dislodged and loose inside the bottle. Do not use bottles if there are any cracks or chips on the bottle. Sudden temperature changes can cause glass bottles to crack or break. Be sure the bottle and liquid are the same temperature when filling or warming. Do not allow baby to self-feed.


Assembly and feeding instructions: Please refer to diagram. Avoid formula bubbles: Do not shake the bottle. Formula should be stirred gently to fully dissolve lumps – as lumps can clog nipple/teat holes. Fill to the desired amount. To prevent potential leaks, do not fill the bottle above the FILL LINE WARNING.

To assemble the internal vent system, connect the reservoir tube to the insert so that they fit snugly together. Place the assembled vent system into the bottle. Insert Infant-Paced Feeding Valve (blue one-way valve) into the nipple/teat. Make sure the valve is fully secured and flush with the nipple/teat base. Insert the nipple/teat into the nipple collar. Affix collar to bottle and turn until closed (do not over-tighten). If warming the bottle – always remove all bottle components before warming and use the storage/travel cap. Do not warm the liquid beyond body temperature of 98.6°F (37°C).

REMINDERS: For proper bottle and venting function, only use Dr. Brown's® nipples/teats with this bottle; other brands may leak or flow too fast or too slow.

When feeding your baby, hold the bottle so that the reservoir tube hole in the bottom of the bottle is not covered by liquid (approximately at a 45° angle). Keep your baby propped up rather than lying down during feeding.

When transporting bottles, always keep them upright and use travel cap or disk to help prevent leaking. During air travel, please remove internal parts and secure with travel cap or disk to prevent leaking from air pressure.

Bottle Components

- 1 Cover/Cap
- 2 Nipple/Teat Collar
- 3 Nipple/Teat
- 4 Infant-Paced Feeding Valve (blue one-way valve)
- 5 Travel Disk
- 6 Insert
- 7 Reservoir Tube
- 8 Vent Assembly: Press firmly together to prevent leakage

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